

Test Sieve Shaker HAVER EML 200 digital plus T with self-readjusting amplitude; G-Control

- For test sieve diameters:
3" (76 mm), 100 mm, 150 mm,
200 mm, 8" (203 mm)
- Programmable, storing functions
- Three-dimensional sieving action
- All electronic control
- Digital display for all functions
- Max. weight of full sieve tower: 3 kg

TECHNICAL DATA:

Operating voltage: 230 V, 50 or 60 Hz,

110 V, 50 or 60 Hz

Amplitude: self-readjusting, max.
apx. 3 mm

Control: programmable

Weight: apx. 34 kg

The EML 200 standard model can accommodate up to 15 half-height test sieves with 200 mm (8", 203 mm) diameter. Sound absorbing cupboard optional.

This test sieve shaker is also available without self-readjusting amplitude in the models HAVER EML 200 digital T and N.



HAVER EML 200 digital plus N
for wet screening with separate control gear