

Non-Electric Models

Used over any effective heat source, the non-electric ALL-AMERICAN® sterilizers provide dry sterile dressings and instruments, with all bacteria and micro-organisms destroyed, in a matter of minutes. Only a small amount of water is needed to create the dry steam.

- UL approved
- One year warranty
- Units include aluminum inner container with rack
- Made in U.S.A.

1915X Features:

- Use over any heat source
- Our most economical unit
- 15 qt/14 liter capacity

1925X Features:

- Use over any heat source
- 46% greater capacity compared to 1915X
- 25 qt/24 liter capacity

1941X Features:

- Use over any heat source
- Our largest unit
- 41 qt/39 liter capacity

SPECIFICATIONS



Model 1915X (15 qt/14 liter)

Gross Capacity	15 qt/14 liter
Overall Height	12 ¹ / ₄ " / 31.2cm
Bottom Height	7 ³ / ₄ " / 19.7cm
Inside Diameter	12 ⁵ / ₈ " / 32.1cm
Unit Weight	15 lbs. / 6.8 kg.
Inner Container No. 2163	
Inside Depth	5 ³ / ₄ " / 14.6cm
Inside Diameter	11 ¹ / ₈ " / 28.3cm
Circumference	35 ⁷ / ₈ " / 91.1cm
Capacity	9.5 qt/9 liter
Volume	550in ³ /9029cm ³
Carton Dimensions	15 ¹ / ₂ " x 14 ¹ / ₂ " x 13 ¹ / ₂ " 39.4cm x 36.8cm x 34.3cm
Shipping Weight	20 lb. / 9.1kg.
Unit Pack: 1	Cube: 1.76

Model 1925X (25 qt/24 liter)

Gross Capacity	25 qt/24 liter
Overall Height	16 ³ / ₄ " / 42.5cm
Bottom Height	12 ¹ / ₄ " / 31.2cm
Inside Diameter	12 ⁵ / ₈ " / 32.1cm
Unit Weight	18 ¹ / ₄ lbs. / 8.3kg
Inner Container No. 2162	
Inside Depth	10 ¹ / ₄ " / 26cm
Inside Diameter	11 ¹ / ₈ " / 28.3cm
Circumference	35 ⁷ / ₈ " / 91.1cm
Capacity	17.6 qt/16.6 liter
Volume	1016in ³ /16.655cm ³
Carton Dimensions	15 ¹ / ₂ " x 14 ¹ / ₂ " x 18 ¹ / ₂ " 39.4cm x 36.8cm x 47cm
Shipping Weight	24 lb. / 10.9kg.
Unit Pack: 1	Cube: 2.41

Model 1941X (41 qt/39 liter)

Gross Capacity	41 qt/39 liter
Overall Height	19" / 48.3cm
Bottom Height	14 ¹ / ₄ " / 36.2cm
Inside Diameter	15 ¹ / ₄ " / 38.7cm
Unit Weight	33 lbs. / 15kg.
Inner Container No. 2164	
Inside Depth	10 ¹ / ₂ " / 26.7cm
Inside Diameter	14" / 35.6cm
Circumference	44 ⁷ / ₈ " / 113cm
Capacity	27.9 qt/26.4 liter
Volume	1613in ³ /26.451cm ³
Carton Dimensions	19" x 19" x 20 ¹ / ₂ " 48.3cm x 48.3cm x 52.1cm
Shipping Weight	41 lb. / 18.6kg.
Unit Pack: 1	Cube: 4.28