

Vortex Mixer

Features

- Suitable for the experiments of mixture, reaction and dissolution of fluid samples in a moment.
- Proper speed according to the environment of experiments can be achieved with a toggle switch.
- The most suitable mixing condition can be made by using a powerful motor.



250VM



260VM

Specifications

Model	250VM	260VM
Type	Test tube	Flat head
Material	NBR	PE foam
Speed	0 - 3000rpm	0 - 3000rpm
Dimension	135W × 185D × 180H(mm)	135W × 185D × 190H(mm)
Net Weight	3 kgs	3 kgs
Power	AC110V or 220V	