

POWERSONIC 400 Series

The Most Powerful Bench-top Ultrasonic Cleaner

POWERSONIC 410 (10liter)



Keypad

■ Features

- Built-in digital thermometer and timer to set optimum condition for cleaning
- Selectable use of ultrasonic power (3 steps : Low-1, Medium-2, High-3)
- Easy to read operation status on digital display
- Built-in drain valve unit
- Silicon rubber heater for durability and stability
- Built-in alarm function for operation stop
- Stainless steel tray (Optional)

■ Specifications



Model	POWERSONIC 405	POWERSONIC 410	POWERSONIC 420
Transducer	Industrial BLT Type		
Time Range	1~99 min (Digital)		
Temp. Range	Ambient~70°C (Digital)		
Frequency	40 KHz		
Housing Material	ABS		
Bath Capacity	5.7 L	10.7 L	20.7 L
Bath Dimensions	W300 x D155 x H150mm	W300 x D240 x H150mm	W500 x D300 x H150mm
Overall Dimensions	W500 x D285 x H260mm	W510 x D375 x H280mm	W720 x D450 x H280mm
Tray Dimensions	W240 x D140 x H150mm	W255 x D200 x H100mm	W500 x D300 x H150mm
Required Power	350 watt	400 watt	700 watt
Net Weight	7 kg	9 kg	14 kg
Voltage	Available from 100 to 240V, 50/60Hz		