

Specifications and designs are subject to change without prior notice for continuous product improvement.

CE RoHS

# POWERSONIC 500 Series

The Most Powerful Bench-top Ultrasonic Cleaner

**HS**  
HWASHIN  
TECHNOLO  
COMPANY

POWERSONIC 510 (10liter)



Keypad

All information available on the "One digital LCD screen"

## ■ Features

- Backlit LCD screen (Furnishing all information)
- Selectable use of ultrasonic power (3 steps : Low-1, Medium-2, High-3)
- Easy to adjust time and temperature setting with number key
- Strong power and durability enhanced by built-in BLT transducer
- Manual function key for turning ON/OFF the heater
- Built-in alarm function for operation stop
- Stainless steel tray (Optional)

## ■ Specifications

Model	POWERSONIC 505	POWERSONIC 510	POWERSONIC 520
Display	70 x 40mm LCD Screen		
Transducer	Industrial BLT Type		
Time Range	1-99 min (Digital)		
Temp. Range	Ambient~70 °C (Digital)		
Frequency	40 KHz		
Housing Material	ABS		
Bath Capacity	5.7 L	10.7 L	20.7 L
Bath Dimensions	W300 x D155 x H150mm	W300 x D240 x H150mm	W500 x D300 x H150mm
Overall Dimensions	W500 x D285 x H260mm	W510 x D375 x H280mm	W720 x D450 x H280mm
Tray Dimensions	W240 x D140 x H150mm	W255 x D200 x H100mm	W500 x D300 x H150mm
Required Power	350 watt	400 watt	700 watt
Net Weight	7 kg	9 kg	14 kg
Voltage	Available from 100 to 240V, 50/60Hz		