



STOMACHER® *accessories & consumables*

➤ Accessories


from Seward

➤ The choice for

Professionals

For maximum effectiveness only genuine Seward Stomacher® bags should be used with your Stomacher® Lab system. Stomacher® bags are extruded from high grade food quality polythene and are guaranteed against bursting. The outstanding strength of the Stomacher® bag ensures safe, quick and efficient handling. The range includes bags for all sizes and types of Stomacher® machine. Filter bags are for processing gritty samples and prevent pipette blocking, closure bags for sample collection and pre-enrichment of processed specimens

complete the bag range. In addition Seward supply Simplettes and bag racks for all bag sizes. Simplettes are convenient, inexpensive disposable pipettes designed to reach the bottom of the Stomacher® 400 bag. They provide a considerable saving over graduated pipettes and pipette tips. All Seward consumable products are supplied irradiated for sterility, certificates of irradiation can be provided for each batch. Only genuine Seward consumables carry the Seward quality logo.





technical

 **Quality Guaranteed**

The unique Seward promise

If any Seward branded Stomacher® bag bursts along the seam during normal laboratory usage, we will supply a replacement box of the same bag free of charge



Product	Model	Type	Code	Bag Dimensions	Quantity	
Bags	Stomacher® 80 Biomaster	Standard Bags	BA6040	152x101mmx65mu	2000	
		Standard Bags	BA6041	177x304mmx65mu	500	
	Stomacher® 400 Classic	Standard Bags Sachets of 5	BA6041/5		500	
		Closure Bags	BA6041/CLR		250	
	Stomacher® 400 Circulator	Strainer Bags	BA6041/STR		200	
		Standard Bags	BA6141	177x304mmx65mu	500	
		Standard Bags Sachets of 5	BA6141/5		500	
		Closure Bags	BA6141/CLR		250	
	Stomacher® 3500	Stomacher® 3500	Strainer Bags	BA6141/STR		200
			Standard Bags	BA6042	380x508mmx65mu	250
Disposable Polypropylene pipettes			BA6051		2000	
Pipettes	Simplettes					
	Bag Racks/ Holders	Stomacher® 80 Biomaster	10 Place Bag Rack	BA6090	1	
		Stomacher® 400 Series	1 Place Bags Rack	BA6091	1	
			10 Place Bag Rack	BA6096	1	
	Stomacher® 3500	1 Place Bag Opener	BA6095	1		
	Stomacher® 3500	3 Place Bag Rack	BA6094	1		

Distributed by:

Seward Limited
98 Great North Road
London N2 0GN
United Kingdom

Telephone: +44 (0) 208 365 4105
Facsimile: +44 (0) 208 365 3999
E-mail: laboratory@seward.co.uk
Web site: www.seward.co.uk

