


Density hydrometers



Ref. temperature: 20 °C

Scale division: 0.002 g/cm³

Error limit: ±0.002 g/cm³

Without thermometer, total length (L): 280 mm

Each packed in a square plastic sleeve

Measuring range g/cm ³	Art.-No.
0.600-0.800	3071001
0.800-1.000	3071002
1.000-1.200	3071003 ^κ
1.200-1.400	3071004 ^κ
1.400-1.600	3071005 ^κ
1.600-1.800	3071006
1.800-2.000	3071007
Complete set with 1 case	3071009
Case for 7 pcs.	201091

^κ Also available as plastic version
(see catalogue page 50)

Available on request with mercury-free
thermometer: 0+35 °C in 1 °C, filling: blue
Error limit: ± 2 °C, total length: 330 mm