

# Lab Burner



Dragon 200

## Specification

Model	Dragon 200
PERFORMANCE DATA	
Max. temperature	1300 °C
Fuel tank capacity	35 g
Working time	60 min
Safe switch	-
Net weight	410 g
Dimension (LxWxH)	8 x 12 x 9 cm

## Features

- ▶ **Reliable & professional**
  - (1) Metal-made base ensures steadiness and spill-proof.
  - (2) Compliant to CE safety regulation for gas-fired appliance.
- ▶ **Stable flame**
  - (1) Superior gas inlet design ensures sufficient butane feed and stable flame outputs.
  - (2) Burns up to 1300°C for max. 60 minutes.
- ▶ **Easy yet safe**
  - Triggered by Piezo Igniter, providing instant flame at one press.

## Ordering Information

**177220-00**

Dragon 220, Lab Burner