

# GL45 Versatile Cap

## Applications



### Waste Bottle

Suitable for collecting waste media by using with GL 45 bottle in connection with suction pump and accessories.




### Moisture Trap Bottle

Available to absorb moisture with desiccant which turns pink when saturated.




### Vapor Absorption (Neutralization) Bottle

Ideal for absorbing acid or alkaline vapor to prevent harmful vapors from spreading into the environment or pump.



### Transfer bottle

Ideal for transferring sterile media stored in GL 45 bottle to another container or reactor with peristaltic pump in a closed system.



### Safety Bottle

Works as a safety bottle to prevent liquid from being sucked into pump.



### Ordering Information

197000-60	Suction Cap for Storage Bottle (1/pk)
197000-60x10	Suction Cap for Storage Bottle (10/pk)
199100-68	Silicone Tube ( $\Phi$ 4.2x $\Phi$ 8 mm), 200cm
197000-60-500	500 ml Storage Bottle with suction cap
197000-60-1000	1000 ml Storage Bottle with suction cap