

Electronic Temperature Regulation Heating Mantle

Summary

It applies to colleges, petrochemical industries, medicine, and environment protection institutions for liquid heating.



Features

- Cold rolling sheet stretching and spraying exterior.
- Electronic temperature regulation mode with characteristics of stepless speed adjustment, uniform temperature, fast heating-up, and safety and durability.
- Featured by anti-corrosion, age-resistant, durable and safe.
- Beautifully shaped, easy to operate and effective.

Specifications

Model	98- I -B									
Specification(ml)	50	100	250	500	1000	2000	3000	5000	10000	20000
Voltage	220V/50Hz									
Max. temperature	0~380°C									
Heating Power (w)	80	100	150	250	350	450	600	800	1200	2500
Working Time	continuous									
Exterior Dimensions W×D×H(mm)	φ 220×165				φ 280×220	φ 330×230	φ 340×245	φ 350×250	φ 425×320	550×510×390
Packaging Dimensions W×D×H(mm)	230×210×185				280×280×275	330×330×280	360×360×280	390×390×300	470×470×380	540×540×420
Net Weight(kg)	2				3.5	4	6	7	10	21
Gross Weight(kg)	2.5				4	5	7	8	12.5	26

Electronic Temperature-regulation Magnetic Stirring Heating Mantle

Summary

It applies to colleges, petrochemical industries, medicine, and environment protection institutions for liquid heating and stirring.



Features

- Cold rolling sheet stretching and spraying exterior.
- Electronic temperature regulation heating mode and magnetic stirring; with characteristics of large heating area, fast warm-up, high temperature evenness, accurate tem-control, time and cost saving.
- Heating and stirring can proceed simultaneously; stepless speed adjustment.

Specifications

Model	98- I -C									
Specification(ml)	50	100	250	500	1000	2000	3000	5000	10000	20000
Voltage	220V/50Hz									
Max. temperature	0~380°C									
Heating Power (w)	80	100	150	250	350	450	600	800	1200	2500
Working Time	continuous									
Exterior Dimensions W×D×H(mm)	φ 220×165				φ 380×220	φ 330×230	φ 340×245	φ 350×250	φ 425×320	550×510×390
Packaging Dimensions W×D×H(mm)	230×210×195				280×280×275	330×330×280	360×360×280	390×390×310	470×470×380	540×540×420
Net Weight(kg)	2.5				5.5	6.5	7.5	8.5	9.8	21
Gross Weight(kg)	2.8				6	7	8	9.2	12	26