

# IKA® Mixing

## Orbital shakers



### KS 130 basic

Small, quiet shaker ensures long life with ideal swivel motion, for a maximum shaking weight of 2 kg.

- Electronic adjustment of speed and timer
- LED display for speed and time adjustment
- Wide range of attachment combinations makes it possible to use almost all shapes and sizes of vessels
- Attachments are not included, please order separately

### Accessories (Page):

Attachments (51): AS 130.1, AS 130.2, AS 130.3, AS 130.4

Shaking motion	orbital
Orbital diameter	4 mm
Max. shaking weight (with attachment)	2 kg
Motor rating input / output	45 / 10 W
Max. ON-time	100 %
Adjustable speed range	80 - 800 rpm
Speed display	LED line
Timer	∞ / 5 - 50 min
Dimensions (W x D x H)	270 x 316 x 98 mm
Weight (without attachment)	7,7 kg
Permissible ambient temperature	5 - 50 °C
Permissible relative humidity	80 %
Protection class acc. to DIN EN 60529	IP 21

Ident. No.  
**2980000** 230 V 50/60 Hz

Ident. No.  
**2980001** 115 V 50/60 Hz