



Kjeldahl Digestion Mantles used for estimation of nitrogen content in fertilisers, dyestuff, drugs etc., provides higher power than the Extraction Mantles.

Kjeldahl Digestion Mantles

The Kjeldahl Digestion Mantle is for use with round bottom or pearshaped Kjeldahl flasks in sizes 15-50 ml.

All Kjeldahl Digestion Mantles will be supplied with a fume tube. Mantles of higher capacity can be supplied on request.

Cat. No.	Cap. ml	No. of flask recess	Total Watts @ 230 V
AV-MK-01	15 - 50	6	720