

Vacuum Filtration System

Two-Year Warranty

Features

- **Special design for large number samples**

MultiVac 310-MS-T is an economical model for large sample filtration. This system includes vacuum pump, manifold, waste bottle and filter holders.
- **Light and handy aluminum design**

The manifold is made of anodized aluminum which is autoclavable, light and cost-effective.
- **Magnetic filter holder is included**

Magnetic filter holder allows for one-handed operation, and can prevent the membrane from being twisted or torn.
- **4 liter PP waste bottle with drain**

The large capacity of 4000ml waste bottle with drain offers convenience of collecting and draining waste solution.



Applications

- Filtration for suspended solid test
- Filtration for microbiology test
- General laboratory filtration

Ordering Information

180311-11(22)-T

MultiVac 310-MS-T, Vacuum Filtration System, complete set includes :

- * Rocker 400, Oil Free Vacuum Pump
- * MultiVac 310-MS, 3 Places Aluminum Manifold
- * MF 3, 300ml Magnetic Filter Holder with lid kid (3 pcs)
- * 4000 ml PP Waste Bottle with drain
- * Silicone Tube (2m)

Specification

Model	MultiVac 310-MS-T	
Voltage	110V/60Hz	220V/50Hz
Max. flow rate	37 l/min	34 l/min
Max. vacuum	- 680 mmHg	
Waste bottle	4000 ml, made of PP with drain	
Number of branch	3	
Filter diameter	47mm	
Funnel capacity	300ml	
Effective filtration area	13.1 cm ²	
Hose barb	5/16 inch (8 mm)	
Manifold Part	Material	
Body	aluminum	
Control valve	aluminum	
Handle	aluminum	
Hose barb	aluminum	