

47 mm Glass Filtration Sets



VF12

Features


- Unique filtration adaptor design can connect with GL45 storage bottle.
- VF12 can filtrate and store solution directly in GL45 storage bottle without the need of being transferred from the filtering flask.

Applications

Ideal for the purification of organic, corrosive liquid such as solvent (mobile phase) for HPLC analysis.

Specifications

Material	VF12
Filter funnel	Borosilicate glass
Funnel support base	Borosilicate glass
Membrane support	Sintered glass
Clamp	Anodized aluminum
Filtration Adaptor	PP (Autoclavable)
Storage Bottle	Borosilicate glass



VF12

Ordering Information

VF12, 47 mm Glass Filtration Set	167200-12	
300 ml Glass filter holder (1+2+3)	167120-31	
1	Glass funnel, 300 ml	167210-03
2	Anodized aluminum clamp	167240-01
3	Glass support base (Including No.8 silicone stopper)	167230-03
4	Filtration Adaptor	197000-65
5	1000ml Storage Bottle	197000-59