

47 mm Stainless Steel Filtration Sets



VF2



VF11

Applications

For filtration requiring the sterilization such as microbiology test or use for purification of corrosive liquid, solvent etc.

Specifications

Material	VF2	VF11
Filter funnel	SUS316	SUS316
Funnel support base	SUS316	SUS316
Membrane support	SUS316	SUS316
Stopper	Silicone rubber (No.16)	Silicone rubber (No.18)
Receiver flask	Borosilicate glass	Borosilicate glass
Drain hose	Rubber (non-autoclavable)	Rubber (non-autoclavable)
Specifications		
Funnel capacity	100 ml	500 ml
Flask capacity	1000 ml	2000 ml
Filter diameter	47 mm / 50 mm	47 mm / 50 mm
Effective filtration area	9.6 cm ²	9.6 cm ²
Hose barb	5/16 inch (8 mm)	5/16 inch (8 mm)

Ordering Information

VF2, 47 mm Stainless Steel Filtration Set 167200-02

100 ml Stainless steel filter holder (1+2+3+4) 167103-10

1	Stainless steel funnel, 100 ml	167103-61
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.16)	167110-17

5	Receiver flask with drain, 1000 ml	167300-08
6	Silicone fixing sucker	167200-45

VF11, 47 mm Stainless Steel Filtration Set 167200-11

500 ml Stainless steel filter holder (1+2+3+4) 167103-50

1	Stainless steel funnel, 500 ml	167103-65
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.18)	167110-18

5	Receiver flask with drain, 2000 ml	167200-13
---	------------------------------------	-----------

