

47 mm Glass Filtration Sets



VF6



VF7

Applications


Ideal for solvent purification and general filtration for microbiology, suspended solid analysis etc.

Specifications

Material	VF6	VF7
Filter funnel	Borosilicate glass	Borosilicate glass
Funnel support base	Borosilicate glass	Borosilicate glass
Membrane support	Sintered glass	S.S. screen with PTFE seal
Clamp	Anodized aluminum	Anodized aluminum
Stopper	Silicone rubber (No.8)	Silicone rubber (No.8)
Receiver flask	Borosilicate glass	Borosilicate glass

Specifications

Funnel capacity	300 ml	300 ml
Flask capacity	1000 ml	1000 ml
Filter diameter	47 mm	47 mm
Effective filtration area	9.6 cm ²	9.6 cm ²
Hose barb	5/16 inch (8 mm)	5/16 inch (8 mm)



VF6

VF7

Ordering Information

VF6, 47 mm Glass Filtration Set		167200-06
300 ml Glass filter holder (1+2+3)		167120-31
1	Glass funnel, 15 ml	167210-03
2	Anodized aluminum clamp	167240-01
3	Glass support base (Including No.8 silicone stopper)	167230-03
4	Receiver flask, 1000 ml	167250-03
5	Silicone fixing sucker	167200-45
VF7, 47 mm Glass Filtration Set		167200-07
300 ml Glass filter holder (1+2+3+4+5)		167120-33
1	Glass funnel, 300 ml	167210-03
2	Anodized aluminum clamp	167240-01
3	Stainless steel support screen	167230-31
4	PTFE gasket	167230-41
5	Glass support base (Including No.8 silicone stopper)	167230-07
6	Receiver flask, 1000 ml	167250-03
7	Silicone fixing sucker	167200-45